



Recombinant Chicken mTERF domain-containing protein 1, mitochondrial (MTERFD1)

Product Code	CSB-EP713721CH-B
Abbreviation	MTERFD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZJC8
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	KGEACA PRSSLPVAVK QEQAAPETVS DLSAAKLCED WDDVPPSSAL EEISEEEAVK IADPLLPPQ SSTLRDYVDH SETLTKLVHL GVDLSQVEKR QKAGQLLLLTL DFEKDITKIL LFLKDVGIED NQLGPFLTKN PYILGEELEA LETRVAYLKS KKFGNAEITQ MVSRAPIYLL FSVERLDNRL GFFKNELGLS VKKTKDLVIR FPRLLTGKLE PVKENLQVCQ VEFGERNEV QQIAFKTPKI LTASKKRLRQ TFDYLNIMG IPHNMLTRFP QVFNSKLLRI KERHMFLIFL GRAQYDPTKP SYISLDQLVS LPDEVFCTEI AKASMQDFEK FLKTI
Source	E.coli
Target Names	MTERF3
Protein Names	Recommended name: mTERF domain-containing protein 1, mitochondrial
Expression Region	65-405
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.