



Recombinant Chicken Transcriptional adapter 2-alpha (TADA2A)

Product Code	CSB-MP713722CH
Abbreviation	TADA2A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZJF3
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MERLASFGND PFDKPPCRGC SSYLTEPYVK CAECGPPPFLL LCLQCFTTRGF EYKKHQSDHT YEIMTSDFPV LDPNWTAEQEE MALLEAVMDC GFGNWQDVAN QMCTKSKEEC EKHYMKHFIN NPLFASTLLN LKQAEAAQHN ETAIPFHPAD DPPRPTFDSL LSRDMAGYMP ARARLSSRSF DNYAEWDLRD IDFVEDDSI LHALKIAVVD IYHSRFKREK TAGRKKIIRD HGLINLRKFQ ILERRYKPEV QDLYETMRRF ARILGPVEHD KFIESHALEF ELRREIKRLQ EYRAAGITNF CSARTYDHLK KTRDEERLKR TMLSEVLQYI QDSSACQQWL SRQADIDSGP TPAAPIPSNS GRRSAPPLNL TGLPGTEKLN EKEKELCQMV RLVPGAYLEY KAALVNECNK QGGLRLAQAR ALIKIDVNKT RKIYDFLIRE GYITKA
Source	Mammalian cell
Target Names	TADA2A
Protein Names	Recommended name: Transcriptional adapter 2-alpha Alternative name(s): Transcriptional adapter 2-like Short name= ADA2-like protein
Expression Region	1-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.