



Recombinant Chicken Protein FAM122A (FAM122A)

Product Code	CSB-MP713759CH
Abbreviation	FAM122A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZLN7
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MAQEKMELDL ELPPGSAAAP SDGGGLRRSN SAPLIHGLSD NSQVFQGSVL RTRRNSTTVM NRHSLFVPSS PIRIPSSRLH QIKQEEGMNL MNRETVHERE VQVAMQMSQS WEESLSLSNN DFEKSSSPKQ VDFVPVSPAP SPTRGIGKQC FPSLQSLVS SSSLPPSPSP SPTRRFSSRR SQSPINCIRP SVLGPIKRKG VTEMEDHPKR LFQGSTNMLS SEALHQPD LG GCLSAHTLDD NRSSAGSSCD SPAEGASTG SPVSLSDSRS PFLPVDLTAK LPID
Source	Mammalian cell
Target Names	FAM122A
Protein Names	Recommended name: Protein FAM122A
Expression Region	1-284
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.