



Recombinant *Methylococcus capsulatus* GTPase Der (der)

Product Code	CSB-BP713860MFM
Abbreviation	der
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q603B5
Product Type	Recombinant Protein
Immunogen Species	<i>Methylococcus capsulatus</i> (strain ATCC 33009 / NCIMB 11132 / Bath)
Purity	≥85% (SDS-PAGE)
Sequence	MLPVVALIGR TNVGKSTLFN YLTRTRDALV ADFPGLTRDR QYGRVQRGER DYFVVDGTTG IEEAEGIDDQ AMKQVDHVLV EADVILFLVD VHAGMTAGDE LIADRLRRIN KPVLLVANKI DGTQALIAGA EFHAFGLGEP LLISAAHGTG VHALLERVEA LLPESGETGE TEAEGGIRIA VVGRPNVGKS TLVNRILGEE RVVVDQPGT TRDSIYIPFE RRGERYTLID TAGIRRRARV NEGVEKFSVI KSFQAIEKAH VVIYLVDAE GLTDQDANLL GMVLEIGRGL LIGFNKWDGL EPEQREKVKR QIDVKLPFLE FAKKYFISAL HGTGVGVLMV AVKPIYQSAM LDLSASRLTQ VLQDCLTAHP PPLVRGRRIK LKYAHQGGHN PPVVVIHGNQ TEDLPAAYRR YLGNEFRNAF KLQGVPLKLV FKSAENPFQG RRNELTERQI RKRRRMIRHV KKR
Source	Baculovirus
Target Names	der
Protein Names	Recommended name: GTPase Der Alternative name(s): GTP-binding protein EngA
Expression Region	1-463
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.