



Recombinant Mouse X-linked lymphocyte-regulated protein 3A (Xlr3a)

Product Code	CSB-MP713873MO
Abbreviation	Xlr3a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60595
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MSSRERKATD TAGRHSRMDP NLSSDDSQNP GAVAAANREV LDAGREDIIS SGTERQQARK EKQDLVQEFE EPPNKVLQEN REKFSRIMTS SFSAMEVKIK DVLKTQCEQR QKLCQDYSLQ FTNLSRKLTS DAYKLKKQAE TLSNMFMEQQ KFIHETLTQ KNRMEEFKSL CEKYLEKLEV LRDSRGNSIA EELRRLIATL EIKLLMLHNQ QNTAAPPQSL LDVLFSS
Source	Mammalian cell
Target Names	Xlr3a
Protein Names	Recommended name: X-linked lymphocyte-regulated protein 3A
Expression Region	1-226
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.