



Recombinant Mouse NF-kappa-B inhibitor beta (Nfkbib)

Product Code	CSB-EP713892MO-B
Abbreviation	Nfkbib
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60778
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MAGVACLGKT ADADEWCDSG LGSLGPDAAA PGGPGLGAEL GPELSWAPLV FGYVTEDGDT ALHLAVIHQH EPFLDFLLGF SAGTEYLDLQ NDLGQTALHL AAILGEASTV EKLYAAGAGV LVAERGGHTA LHLACRVRAH TCACVLLQPR PSHPRDASDT YLTQSQDCTP DTSHAPAAVD SQPNPENEEE PRDEDWRLQL EAENYDGHTP LHVAVIHKDA EMVRLLRDAG ADLNKPEPTC GRTLPLHLAVE AQAASVLELL LKAGADPTAR MYGGRTPLGS ALLRPNPILA RLLRAHGAPPEDEDDKLSP CSSSGSDSDS DNRDEGDEYD DIVVHSGRSQ NRQPPSPASK PLPDDPNPA
Source	E.coli
Target Names	Nfkbib
Protein Names	Recommended name: NF-kappa-B inhibitor beta Short name= NF-kappa-BIB Alternative name(s): I-kappa-B-beta Short name= IκB-B Short name= IκB-beta Short name= IκappaBbeta
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.