



Recombinant Mouse Krueppel-like factor 2 (Klf2)

Product Code	CSB-MP713898MO
Abbreviation	Klf2
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q60843
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MALSEPI LPS FATFASPCER GLQERWPRNE PEAGGTDEDL NNVLDFILSM GLDGLGAENP PEP PPQPPPP AFYYPEPGAP PPYSIPAASL GTELLRPDLD PPQGPALHGR FLLAPPGLRV KAEPPEVDGG GYGCAPGLAH GPRGLKLEGA PGATGACMRG PAGRPPPPPD TPPLSPDGPL RIPASGPRNP FPPPFPGGPS FGGPGPALHY GPPAPGAFGL FEDAAAAAAA LGLAPPATRG LLTPPSSPLE LLEAKPKRGR RSWPRKRAAT HTCSYTNC GK TYTKSSHLKA HLRHTHTGEKP YHCNWE GCGW KFARSDELTR HYRKHTGHRP FQCHLCDRAF SRSDHLALHM KRHM</p>
Source	Mammalian cell
Target Names	Klf2
Protein Names	Recommended name: Krueppel-like factor 2 Alternative name(s): Lung krueppel-like factor
Expression Region	1-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>