



Recombinant Spindle assembly abnormal protein 6 (sas-6)

Product Code	CSB-MP713948CXX
Abbreviation	sas-6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60P76
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis briggsae
Purity	≥85% (SDS-PAGE)
Sequence	<p> MTSRLALFDQ TLSASLLQRV NSDPHDYKAY KTRVCLKILE QRGEHGKEL KFEISRADDF EFLFAAILNN EKYQELAKAH ELTVDFDFTP KVIIQHLLCK NIVKTADEDG IEVDARAKPG YHSIIENRGV PPTEINIKLE KDKNSCEFEI FSKTPISKGK IFAMSLQAVR GDLLIAHLLK ICSSQSTKIA TLHRSTDELA ALKTKVAEME KRNQSIEDLE KKIEEQEERI REMEDELELV KEERENIRVI AEEKEDLVAE LQNDLESVNR ELEENQEELE IVGKMLGEEQ GKVDQLQKRN VAHQKELATL QRNFQKVESQ LKRSSVEQSQ QSLDIRKLRE LEADLKEKDS MVENLTQTVG ILKKELEDER VKMQETMTSF ESLKVENESV KERLAMYRNQ RFSPAPVGLQ TTTIPSAGVG FKTVLGQNSP YANNPNMRTP FRDTSSALNT GLTYQNSITP LPNTRFLMA DDTTGSSVTN TPPVLRNPL </p>
Source	Mammalian cell
Target Names	sas-6
Protein Names	Recommended name: Spindle assembly abnormal protein 6
Expression Region	1-489
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.