



Recombinant Putative elongator complex protein 4 (elpc-4)

Product Code	CSB-BP713954CXX
Abbreviation	elpc-4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60ZM2
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis briggsae
Purity	≥85% (SDS-PAGE)
Sequence	MLNIGDSVQI PGCSTKKRLL ETTSGCDSFD TLLGGALVNS SIVLIDEYRS RCYGSYLVR S FLAEG LHYGH RCFVADPSED VSETILKVIP TRKAAEDDQK KEIPSSPGQA DNDMKIAWRY GNVKQVTSEI GPSGNQYDFS KHLENSNVDV YNEEIPSFSG LYKILCEVVQ KEEAHTKSTG RGGPKKNLLR VVLKNIHLNI WEDYESLGRF LCCLRSLARS SYAVIYIVAN SHRMPNDTWR QLESAADTHI QLMPFNEQEK KMFRHLGTAH GYFHLKSLPR LMSAGTHTPP ILDLIFEASS KKG FQIRVMH LPPAFDEPAP GQQKATACQN IDF
Source	Baculovirus
Target Names	elpc-4
Protein Names	Recommended name: Putative elongator complex protein 4 Short name= ELP4
Expression Region	1-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.