



Recombinant Mouse Transcription factor AP-2 gamma (Tfap2c)

Product Code	CSB-MP713969MO
Abbreviation	Tfap2c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q61312
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MLWKITDNVK YEEDCEDRHD SSSNGNPRIP HLSSPGQHLY SPAPPLSHTG VAEYQPPPYF PPPYQQLAYS QSADHYSHLG EAYAAAMNPL HQPAATGSQQ QAWPGRQSQE GSSLASHHSR SASLIPHISG LEGGSVSARR EVYRRSDLLL PHAHALEAGL AENLGLHEMA HPIEEVQNVD DAHLLLHDQT VIRKGPISMT KNPLGLPCQK DLVGVMNPS EVFCSVPGRL SLLSSTSKYK VTVAEVQRR SPPECLNASL LGGVLRRAKS KNGGRSLREK LDKIGLNLP A GRRKAAHVTL LTSLVEGEAV HLARDFAYVC EAEFPSKAVA DYLTRPHLGG RNEMATRKSM LLAAQQVCKE FTDLLHQDRT PNGNNRPAQV LEPNIQNCLS HFSLITHGFG SQAICAAVSA VQNYIKEALI AIDKSYMNP G DQSPADSSKT MEKMEKHRK
Source	Mammalian cell
Target Names	Tfap2c
Protein Names	Recommended name: Transcription factor AP-2 gamma Short name= AP2-gamma Alternative name(s): AP-2.2 Activating enhancer-binding protein 2 gamma
Expression Region	1-449
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is a sequence-specific DNA-binding transcription factor involved in the activation of several developmental genes. The encoded protein can act as either a homodimer or heterodimer with other family members and is induced during retinoic acid-mediated differentiation. It plays a role in the development of the eyes, face, body wall, limbs, and neural tube.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.