



Recombinant Mouse Double-strand-break repair protein rad21 homolog (Rad21)

Product Code	CSB-EP713978MO-B
Abbreviation	Rad21
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q61550
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MFYAHFVLSK RGPLAKIWLA AHWDKCLTKA HVFECNLESS VESIISPKVK MALRTSGHLL LGVVRIYHRK AKYLLADCNE AFIKIKMAFR PGVVDLPEEN REAYNAITL PEEFHDFDQP LPDLDDIDVA QQFSLNQRSV EEITMREEVG NISILQENDF GDFGMDDREI MREGSAFEDD DMLVSTSASN LLEPEQSTS NLNEKMNHLE YEDQYKDDNF GEGNDGGILD DKLISNNDGG IFDDPPALSE AGVMLPEQPA HDDMDEDDNG SLGGPDSPDS VDPVEPMPTM TDQTTLVPNE EEAFALEPID ITVKETKAKR KRKLIVDSVK ELDSKTIRAQ LSDYSDIVTT LDLAPTKKL MMWKETGGVE KLFFLPAQPL WNNRLLKLFT RCLTPLVPED LRKRRKGGEA DNLDEFLKEF ENPEVPREEQ QPQQQQPQPQ RDVIDEPIIE EPSRLQDSVM EASRTTIEES AMPPPPPQGV KRKAGQIDPE PSIPPQQVEQ MEIPPVELPP EEPNICQLI PELELLPEKE KEKEKEKEEE EEEEDEDASG GDQDQEERRW NKRTQQMLHG LQRALAKTGA ESISLLELCR NTNRKQAAAK FYSFLVLKKQ QAIELTQEEP YSDIIATPGP RFHII
Source	E.coli
Target Names	Rad21
Protein Names	Recommended name: Double-strand-break repair protein rad21 homolog Alternative name(s): Pokeweed agglutinin-binding protein 29 Short name= PW29 SCC1 homolog
Expression Region	1-635
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is highly similar to the gene product of Schizosaccharomyces pombe rad21, a gene involved in the repair of DNA double-strand breaks, as well as in chromatid cohesion during mitosis. This protein is a nuclear phospho-protein, which becomes hyperphosphorylated in cell cycle M phase. The highly



regulated association of this protein with mitotic chromatin specifically at the centromere region suggests its role in sister chromatid cohesion in mitotic cells.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.