



Recombinant Mouse GTP-binding nuclear protein Ran, testis-specific isoform (Rasl2-9)

Product Code	CSB-MP713989MO
Abbreviation	Rasl2-9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q61820
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	AAQGEPQVQ FKVVLVGDGG TGKTTFMKRH LTGEFEKEYV ATLGVEVHTL VFHTNRGPIK FNWDTAGQE KFGGLRDGYY IQAQCAIMF DVTSRVTYKN VPSWHKDLVR VCENIPIVLC GNKVDVKDMK VKAKPILFHR KKNLQYYDIS ARSNYNFEKP FFWLARKLIG DPNLEFVAMP ALAPPEVVMD PALAAQYEHD LEVAQTTALP DEEDDL
Source	Mammalian cell
Target Names	Rasl2-9
Protein Names	Recommended name: GTP-binding nuclear protein Ran, testis-specific isoform
Expression Region	2-216
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.