



# Recombinant Probable V-type proton ATPase subunit B (vha-12)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP713993CXX                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Abbreviation</b>      | vha-12                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                      |
| <b>Uniprot No.</b>       | Q619C0                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Immunogen Species</b> | Caenorhabditis briggsae                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Sequence</b>          | MAAVDVNQTI TGHKSAIIRN YTTNPRLIQ TVSGVNGPLV ILNDVKFPQF<br>SEIVKITLPD GSQRSGQVLE IAKNKAVVQV FEGTSGIDAK NTICEFTGDI<br>LRTPVSEDML GRIFNGSGKP IDKGPPVLA E DFLDINGQPI NPWSRIYPEE<br>MIQTGISAID VMNSIARGQK IPIFSAAGLP HNEIAAQIVR QGGLVQLPDR<br>PHEQTNFAIV FAAMGVNMET ARFFKQDFEE NGSMENVCLF LNLANDPTIE<br>RIITPRIALT AAEFLAYQCS KHLVVLVLTDM SSYAEALREV SAAREEVPGR<br>RGFPGYMYTD LATIYERAGR VEGRDGSITQ IPILTMPNDD ITHPIPDLTG<br>YITEGQIYVD RQLHNRLIYP PINVLPSLSR LMKSAIGDKM TREDHSDVSN<br>QLYACYAIGK DVQAMKAVVG EEALSSDDL YLEFLVKFEK NFISQGHYEN<br>RTIFESLDIG WELLRIFPRE MLKRIPAKSL DKYYPRGGAK E |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Target Names</b>      | vha-12                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Protein Names</b>     | Recommended name: Probable V-type proton ATPase subunit B Short name= V-ATPase subunit B Alternative name(s): Vacuolar proton pump subunit B                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Expression Region</b> | 1-491                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                                             |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life                                                                                                                                                                                                                                                                                                                                                       |



of lyophilized form is 12 months at -20°C/-80°C.