



Recombinant Guanine nucleotide-binding protein subunit beta-1 (gpb-1)

Product Code	CSB-BP714006CXX
Abbreviation	gpb-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q61ZF6
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis briggsae
Purity	>85% (SDS-PAGE)
Sequence	MSELDQLRQE AEQLKSQIRE ARKSANDTTL ATVASNLEPI GRIQMRTRRT LRGHLAKIYA MHWASDSRNL VSASQDGKLI VWDSYTTNKV HAIPLRSSWV MTCAYAPSGS FVACGGLDNI CSIYSLKTRE GNVRVSRELP GHTGYLSCCR FLDDNQIVTS SGMTCALWD IETGQQCTAF TGHTGDVMSL SLSPDFRTFI SGACDASAKL WDIRDGMCKQ TFPGHESDIN AVAFFPSGNA FATGSDDATC RLFDIRADQE LAMYSHDNII CGITSVAFSK SGROLLFAGYD DFNCNVWDSM RQERAGVLAV HDNRVSLGLV TEDGMAVCTG SWDSFLKIWN
Source	Baculovirus
Target Names	gpb-1
Protein Names	Recommended name: Guanine nucleotide-binding protein subunit beta-1
Expression Region	1-340
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.