



Recombinant Mouse Keratin, type I cuticular Ha2 (Krt32)

Product Code	CSB-YP714011MO
Abbreviation	Krt32
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62168
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MPSVCMPTTY RPASCLSKTY LSSSCQPSNR RPTGCISSSM GTYGLFCEGA FNGNEKETMQ VLNDRLANYL EKVRQLEKEN AELEGKIQDV YQGQVLT MCP DYQSYFQTIE ELQQKVLCTK AENARMIVHI DNAKLAADD F RTKYETELAL RQLVEADTNG LRRILDELTL NKADLEAQVE SLKEELLCLK RNHEEEVGV L RQQLGDR LNI EVDAAPPVDL TRMLEEMRCQ YETMVETNHR DVEEWFNMQM EELNKQVATS SEQLQSYQSD IIDLRRTVNT LEIELQAQHS LRDSLENTLG ETEGRFTSQL SQMQCMITNV ESQ LSDIRCD LERQNQEYKV LLDVKARLEC EIDTYRGLLE SEDSKLPCNP CSTPSCQPCA PSPGVSR TVC VPHTVCVPCS PCLQTRY
Source	Yeast
Target Names	Krt32
Protein Names	Recommended name: Keratin, type I cuticular Ha2 Alternative name(s): Hair keratin, type I Ha2 Keratin-32 Short name= K32
Expression Region	1-407
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is a member of the keratin gene family. As a type I hair keratin, it is an acidic protein which heterodimerizes with type II keratins to form hair and nails. The type I hair keratins are clustered in a region of chromosome 17q12-q21 and have the same direction of transcription.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.