



# Recombinant Mouse Neuronal pentraxin-1 (Nptx1)

<b>Product Code</b>	CSB-BP714021MO
<b>Abbreviation</b>	Nptx1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q62443
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	QDFGPTRF ICTSVPVDAD MCAASVAAGG AEELRSNVLQ LRETVLQQKE TILSQKETIR ELTTKLGRCE SQSTLDSPGP EARS GGGRKQ PGSGKNTMGD LSRTPAAETL SQLGQTLQSL KTRLENLEQY SRLNSSSQTN SLKDLLQSKI DDLERQVLSR VNTLEEGKGG PKNDTEERAK IESALTS LHQ RISELEKGGQK DNRPGDKFQL TFPLRTNYMY AKVKKSLPEM YAFTVCMWLK SSAAPGVGTP FSYAVPGQAN ELVLI EWGNN PMEILINDKV AKLPFVINDG KWHHICVTWT TRDGVWEAYQ DGTQGGNGEN L APYHPIKPQ GVLVLGQE QD TLGGGFDATQ AFV GELAHFN IWDRKLTPGE VYNLATCSSK ALSGNVIAWA ESQIEIFGGA TKWTFEACRQ IN
<b>Source</b>	Baculovirus
<b>Target Names</b>	Nptx1
<b>Protein Names</b>	Recommended name: Neuronal pentraxin-1 Short name= NP1 Alternative name(s): Neuronal pentraxin I Short name= NP-I
<b>Expression Region</b>	23-432
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	NPTX1 is a member of the neuronal pentraxin gene family. Neuronal pentraxin 1 is similar to the rat NP1 gene which encodes a binding protein for the snake venom toxin taipoxin. Human NPTX1 mRNA is exclusively localized to the nervous system.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.