



Recombinant Rat Galectin-8 (Lgals8)

Product Code	CSB-BP714024RA
Abbreviation	Lgals8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62665
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	<p> MLSLSNLQNI IYNPTIPYVS TITEQLKPGS LIVIRGHVPK DSERFQVDFQ HGNSLKPRAD VAFHFNPRFK RSNCIVCNTL TNEKWGWEEI THDMPFRKEK SFEIVIMVLK NKFHVAVNGK HILLYAHRIN PEKIDTLGIF GKVNIHSIGF RFSSDLQSME TSTLGLTQIS KENIQKSGKL HLSLPFEARL NASMGPGRTV VVKGEVNTNA TSFNVDLVAG RSRDIALHLN PRLNVKAFVR NSFLQDAWGE EERNITCFPF SSGMYFEMII YCDVREFKVA VNGVHSLEYK HRFKDLSSID TLAVDGDIRL LDVRSW </p>
Source	Baculovirus
Target Names	Lgals8
Protein Names	Recommended name: Galectin-8 Short name= Gal-8 Alternative name(s): 30 kDa S-type lectin RL-30
Expression Region	1-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a member of the galectin family. Galectins are beta-galactoside-binding animal lectins with conserved carbohydrate recognition domains. The galectins have been implicated in many essential functions including development, differentiation, cell-cell adhesion, cell-matrix interaction, growth regulation, apoptosis, and RNA splicing. This gene is widely expressed in tumoral tissues and seems to be involved in integrin-like cell interactions. Alternatively spliced transcript variants encoding different isoforms have been identified.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.