



Recombinant Bacillus cereus Glycogen synthase (glgA)

Product Code	CSB-EP714123BAAD
Abbreviation	glgA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q632H4
Product Type	Recombinant Protein
Immunogen Species	Bacillus cereus (strain ZK / E33L)
Purity	≥85% (SDS-PAGE)
Sequence	MNILFAVSEC VPFVKSGGLA DVAGALPKEL KKLGVVEVRIL LPNYSILIPQK LRDGCTLHKV INVPLGWRNQ YCGILKGEQD GITYYLIDNE YYFKRDSLYG HYDDGERFSY FSKAVLECIP HLD FEVDVLH SHDWHTAMVN FLLREKYQDN PLYEHIKTVY TIHNLQFQGV FPPEVMYDLL ELGDEYFHSE QLEFYGNVNF MKGGIIASDQ ITAVSPTYKE EIQYEFFGEK LDGLLRKYND KLSGIVNGID TSVYNPETDS YITAQYDAGS LYKNNENKRA LQRYFGLPEK EDTPIISMVT RLTKQKGLDL VRTVFREIME EDVQCIILGS GDSEYEQFFE WMAYEYPEKV KVYIGFNEEL AHQVYAGSDL FLMPSLFEPC GLGQLIALAY GTIPIVRETG GLNDTVQSYD EETGEGNGFS FTNFNAHDML HTVLRAIEFY HDKPVWEQLV KQAMTEDYSW EKSALAYKKL YKGLME
Source	E.coli
Target Names	glgA
Protein Names	Recommended name: Glycogen synthase EC= 2.4.1.21 Alternative name(s): Starch [bacterial glycogen] synthase
Expression Region	1-476
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.