



# Recombinant Rat Lipopolysaccharide-binding protein (Lbp)

<b>Product Code</b>	CSB-MP714125RA
<b>Abbreviation</b>	Lbp
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q63313
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	VNPAM VVRITDKGLE YAAKEGLLSL QRELYKITLP DFSGDFKIKAVGRGQYEFHS LEIQSCQLRG SSLKPLPGRG LSLISDSSI SVRGKWKVRR SFVKLHGSFD LDVKSVTISV DLLLGVDPSE RPTVTASGCS NRIRDLELHV SGNVGWLLNL FHNQIESKLQ KVLESKICEM IQKSVTSDLQ PYLQTLPVTA DIDTILGIDY SLVAAPQAKA QTLDMFKGE IFNRNHRSPV TTPTPTMSLP EDSKQMVYFA ISDQAFNIAT RVYHQAGYLN FTITDDMLPP DSNIRLNTKA FRPFTPLITR KYPDMNLELL GTVVSAPLLN VSPGNLSLAP QMEIEGFVIL PSSARESVFR LGVVTNVFVS LTFDNSKVTG MLHPEKAQVR LIESKVGGMFN VNLFQAFLNY YLLNSLYPDV NDELAKGFPL PLPRRIKLHD LDFQIHKNFL YLGANVQYMR V
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Lbp
<b>Protein Names</b>	Recommended name: Lipopolysaccharide-binding protein Short name= LBP
<b>Expression Region</b>	26-481
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein is involved in the acute-phase immunologic response to gram-negative bacterial infections. Gram-negative bacteria contain a glycolipid, lipopolysaccharide (LPS), on their outer cell wall. Together with bactericidal permeability-increasing protein (BPI), the encoded protein binds LPS and interacts with the CD14 receptor, probably playing a role in regulating LPS-dependent monocyte responses. Studies in mice suggest that the encoded protein is necessary for the rapid acute-phase response to LPS but not for the clearance of LPS from circulation. This protein is part of a family of structurally and functionally related proteins, including BPI, plasma cholesteryl ester transfer



protein (CETP), and phospholipid transfer protein (PLTP). Finally, this gene is found on chromosome 20, immediately downstream of the BPI gene.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.