



Recombinant *Xenopus laevis* Protein FAM136A (fam136a)

Product Code	CSB-BP714457XBE
Abbreviation	fam136a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q63ZH8
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	≥85% (SDS-PAGE)
Sequence	MAEEHQNRLQ NAVDSMVKSL ERDNIRKMQG KMFRCSAQCC EDNGASMQQV HHCIERCHTP LAQAQSLVTT ELERFQNRLS RCTMNCSDQA KDAFDSGSKE AQVKAQLEGC VIKCAEEHMN LIPSMTKKLK DALAQADK
Source	Baculovirus
Target Names	fam136a
Protein Names	Recommended name: Protein FAM136A
Expression Region	1-138
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.