



Recombinant Rat EH domain-containing protein 1 (Ehd1)

Product Code	CSB-MP714475RA
Abbreviation	Ehd1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q641Z6
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MFSWVSKDAR RKKEPELFQT VAEGLRQLYA QKLLPLEEHY RFHEFHSPAL EDADFDNKPM VLLVGQYSTG KTFIRHLIE QDFPGMRIGP EPTTDSFIIV MHGPTGVVP GNALVVDPRR PFRKLNAGFN AFLNRFMCAQ LPNPVLDSIS IIDTPGILSG EKQRISRGYD FAVLEWFAE RVDRIILLFD AHKLDISDEF SEVIKALKNH EDKIRVVLNK ADQIETQQLM RYVYALMWSL GKIINTPEVV RVYIGSFWSH PLLIPDNRKL FEAEQDLFK DIQSLPRNAA LRKLNLIKR ARLAKVHAYI ISSLKKEMPN VFGKESKKKE LVNNLGEIYQ KIEREHQISS GDFPSLRKMQ ELLQTQDFSK FQALKPKLLD TVDDMLANDI ARLMVMVRQE ESLMPSQAVK GGAFDGTMNG PFGHGYGEGA GEGIDDVEWV VGKDKPTYDE IFYTLSPVNG KITGANAKKE MVKSKLPNTV LGKIWKLADV DKDGLLDDEE FALANHLIKV KLEGHELPAD LPPHIPPSPK RRHE
Source	Mammalian cell
Target Names	Ehd1
Protein Names	Recommended name: EH domain-containing protein 1
Expression Region	1-534
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene belongs to a highly conserved gene family encoding EPS15 homology (EH) domain-containing proteins. The protein-binding EH domain was first noted in EPS15, a substrate for the epidermal growth factor receptor. The EH domain has been shown to be an important motif in proteins involved in protein-protein interactions and in intracellular sorting. This protein is thought to play a role in the endocytosis of IGF1 receptors.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.