



Recombinant Fowl adenovirus A serotype 1

Protein GAM-1 (8)

Product Code	CSB-MP714520FCL
Abbreviation	8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64770
Product Type	Recombinant Protein
Immunogen Species	Fowl adenovirus A serotype 1 (strain CELO / Phelps) (FAdV-1) (Avian adenovirus gal1 (strain Phelps))
Purity	≥85% (SDS-PAGE)
Sequence	MARNPFRMFP GDLPYMGTI SFTSVVPVDP SQRNPTTSLR EMVTTGLIFN PNLTGEQLRE YSFSPLVSMG RKAIFADYEG PQRIIHVTIR GRSAEPKTPS EALIMMEKAV RGAFAPDWV AREYSDPLPH GITHVGD LGF PIGSVHALKM ALDTLKIHPV RGVGVPGYEG LCGTTTIKAP RQYRLTTGV FTKKDLKRTL PEPFFSRFFN QTPEVCAIKT GKNPFSTEIW CMTLGGDSPA PERNEPRNPH SLQDWARLGV METCLRMSRR GLGSRHHPYH SL
Source	Mammalian cell
Target Names	8
Protein Names	Recommended name: Protein GAM-1 Alternative name(s): Gallus-anti morte protein Gam1
Expression Region	1-282
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.