



Recombinant Apple stem pitting virus Capsid protein (ORF5)

Product Code	CSB-MP714523AEAA
Abbreviation	ORF5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64966
Product Type	Recombinant Protein
Immunogen Species	Apple stem pitting virus (isolate PA66) (ASPV)
Purity	≥85% (SDS-PAGE)
Sequence	MTSNGSQPQA STPMVSAEEP AAAASVPNST PMVSAEGPAA AVSAPNSSVV SSAPASAPTA SEPVISQVQS LAPIVSGFDP NLHGRLTNEQ MRQAQNEAAM QGYEEGSRRN PRLPSSTTAH NDYASMNSNP FETGTAYGGA PRVSFGSYPT FPGSGSASEP NSQRIFPQQH GVNPPAHASD LVPHQATSGG NTGTPFTLGN RAPRNATANT GGMRRRLDSV GLKNIRYEPQ AGVVASNQKI RAVGVALIGM GIPEHQLTEV GVYLARHCAD VGASDKSALL GTFPGSDITL EEVGTMIKQT EGCTLRQYCA FYAKHVWNLM LQTQSPPANW VGKEFKFETR YAAFDFFFGV ESTASLEPAD GLIRLPTQAE RVANATSKEI QMYRIRSMEG TQAVNFG EVT GGKIGPKPVL SIRK
Source	Mammalian cell
Target Names	ORF5
Protein Names	Recommended name: Capsid protein Short name= CP Alternative name(s): Coat protein
Expression Region	1-414
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.