



Recombinant Chicken Transcription factor Sp8 (SP8)

Product Code	CSB-EP714528CH
Abbreviation	SP8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64HY5
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MLAATCNKIG SPSPSPSALS DSASSFGKGF HPWKRSSSSS SASAGSCGAV GSGLPFGVA GAARNGSSAA AAAAAAAAAA LVSDSFSCGG SPGSSAFSLT SSGAAAASSP FANDYSVFQA PGSAGGGGGG GGGGGGAAGQ EAAHQPVFIS KVHASVEGLQ GIYPRVGM AH PYESWFKPSH PGLAAGEVGS AGASSWWDVG AGWIDVQSPN GAAALPGSLH PAAGGLQSSL HSPLGGYNSD YSGLGHSAFG GGASSHLLSP AGQHLMDFGK PVLPGSYPDS APSPLAGAGG SMLGGGPAAP LSASPRSSAR RYSGRATCDC PNCQEAERLG PAGASLRRKG LHSCHIPGCG KVYGKTSHLK AHLRWHTGER PFVCNWLFCG KRFRSDELQ RHLRHTTGEK RFACPVCNKR FMRSDHLSKH VKTHSGPGGA GPGGGGGPGP GPGGKKGSDT DSEHSAAGSP PCHSPELLPP PEPGHRNGLE
Source	E.coli
Target Names	SP8
Protein Names	Recommended name: Transcription factor Sp8
Expression Region	1-480
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.