



Recombinant *Oryza sativa* subsp. *japonica* Probable protein phosphatase 2C 55 (Os06g0526700, LOC_Os06g33530)

Product Code	CSB-BP714591OFG
Abbreviation	Os06g0526700, LOC_Os06g33530
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q652Z7
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MRRHLLGLL RRAAASSTSA ASSRAGPHPS LHAPGPLRNG GSAPRFFSSR GGAGAASKGL GDDEVELYSL LLGVSIGDEG EASSRGPAAAS RGGRRGRNSK RQPPRSRFDG DGVGCSKDGK LSWGYSFQG RRPSMEDRLS IKSTTVNGET VSLFGVFDGH GGPRAAEYLK KHLFKNLVKH PKFLKDTKLA INQTFKTD ADFLQSISSDR YRDDGSTAVA AILIGNRLYV ANVGDSRAVA LKAGKAVPLS EDHKPNKKDE RKRIEDAGGI VVSDDIWRVD GILAVSRAFG NRLMKRYVKA EPNIQEKVVD EGGLEYLVLAT DGLWDVMRNE DAVSLLKAQD GPKAAAMKLT EVARSRLTLD NVTCIVLQFH HGKSTNSK
Source	Baculovirus
Target Names	Os06g0526700
Protein Names	Recommended name: Probable protein phosphatase 2C 55 Short name= OsPP2C55 EC= 3.1.3.16
Expression Region	1-378
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.