



Recombinant *Bacillus licheniformis* Ribose import ATP-binding protein RbsA (rbsA)

Product Code	CSB-MP714611BQU
Abbreviation	rbsA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q65E55
Product Type	Recombinant Protein
Immunogen Species	<i>Bacillus licheniformis</i> (strain ATCC 14580 / DSM 13 / JCM 2505 / NBRC 12200 / NCIMB 9375 / NRRL NRS-1264 / Gibson 46)
Purity	≥85% (SDS-PAGE)
Sequence	MHIKMKDIFK AFGQNQVLSG VSFELQEGEV HALMGENGAG KSTLMNLLTG LHKLDSTGIE IDGKETYFSD PKEAEQNGIA FIHQELNIWP EMTVLENLFI GRELSSKLG F LNNKKMKALA KEQLERLGVS ISLEKEAGDC SVGQQQMIEI AKALMTDAKV IIMDEPTAAL TEREIEKLF G VIRALKKNGV SIVYISHRME EIFTICDRIT VMRDGKTVDT KRIPDTDFHE VVKKMGREL TERYPERQPN PGRVVLEV KQ ASKKG V FQNI SFSVRSGEIV GISGLMGAGR TELMRAVFGL DPLDSGDIFI EGKKAALKKP SDAVQKGIGF ITENRKDEGL VLDTSIRENI ALPNLASFSP KGWIDKKSEQ EFVDLLIKRL TIKTESPETH ARNLSSGGNQQ KVVIKWIGI GPKVLILDEP TRGVDVGA KR EIYQLMNE LT DRGVAIVMVS SELPEILGMS DRVLVIHEGT LSGELSRNDA TQERIMTLAT GGR
Source	Mammalian cell
Target Names	rbsA
Protein Names	Recommended name: Ribose import ATP-binding protein RbsA EC= 3.6.3.17
Expression Region	1-493
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.