



# Recombinant *Bacillus licheniformis* GTPase obg (obg)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP714626BQU                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Abbreviation</b>      | obg                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                |
| <b>Uniprot No.</b>       | Q65GM7                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Immunogen Species</b> | <i>Bacillus licheniformis</i> (strain ATCC 14580 / DSM 13 / JCM 2505 / NBRC 12200 / NCIMB 9375 / NRRL NRS-1264 / Gibson 46)                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Sequence</b>          | MFVDQVKIYV KGGDGGNGMV AFRREKYVPK GGPAGGDGGK<br>GGDVVFKVDE GLSTLMDFRY QRHFKAARGE HGMSKNQHGR<br>NAEDMVVKVP PGTVIDDDT KQVIADLTEH GQEAVIAKGG RGGRGNTRFA<br>TPANPAPQLS ENGEPGKERY IVLELKV LAD VGLVGFP SVG KSTLLSVVSS<br>AKPKIADYHF TTLNPNLGMV ETEDGRSFVM ADLPGLIEGA HEGVGLGHQF<br>LRHIERTVI VHVIDMSGLE GRDPYEDYVT INKELEQYNL RLTERPQIIV<br>ANKMDMPDAE ENLKAFKEKL TDDYPVFPIS AVTRQGLRDL LFEIADRLET<br>TPEFPLYDEE DMAENRV MYK LEDEEAPFEI SRDPDGT FVL SGAKLERLFK<br>MTDFSRDESV KRFARQLRGM GVDDALRARG AKDGD TIRLL EFEFEFID |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Target Names</b>      | obg                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Names</b>     | Recommended name: GTPase obg Alternative name(s): GTP-binding protein obg                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Expression Region</b> | 1-428                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                       |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life                                                                                                                                                                                                                                                                                 |



of lyophilized form is 12 months at -20°C/-80°C.