



Recombinant Bacillus licheniformis GTPase obg (obg)

Product Code	CSB-MP714626BQU
Abbreviation	obg
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q65GM7
Product Type	Recombinant Protein
Immunogen Species	Bacillus licheniformis (strain ATCC 14580 / DSM 13 / JCM 2505 / NBRC 12200 / NCIMB 9375 / NRRL NRS-1264 / Gibson 46)
Purity	≥85% (SDS-PAGE)
Sequence	MFVDQVKIYV KGGDGGNGMV AFRREKYVPK GGPAGGDGGK GGDVVFKVDE GLSTLMDFRY QRHFKAARGE HGMSKNQHGR NAEDMVVKVP PGTVIDDDT KQVIADLTEH GQEAVIAKGG RGGRGNTRFA TPANPAPQLS ENGEPGKERY IVLELKV LAD VGLVGFP SVG KSTLLSVVSS AKPKIADYHF TTLNPNLGMV ETEDGRSFVM ADLPGLIEGA HEGVGLGHQF LRHIERTVI VHVIDMSGLE GRDPYEDYVT INKELEQYNL RLTERPQIIV ANKMDMPDAE ENLKAFKEKL TDDYPVFPIS AVTRQGLRDL LFEIADRLET TPEFPLYDEE DMAENRVMYK LEDEEAPFEI SRDPDGT FVL SGAKLERLFK MTDFSRDESV KRFARQLRGM GVDDALRARG AKDGD TIRLL EFEFEFID
Source	Mammalian cell
Target Names	obg
Protein Names	Recommended name: GTPase obg Alternative name(s): GTP-binding protein obg
Expression Region	1-428
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.