



Recombinant *Bacillus licheniformis* Allantoinase (allB)

Product Code	CSB-MP714658BQU
Abbreviation	allB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q65LN0
Product Type	Recombinant Protein
Immunogen Species	<i>Bacillus licheniformis</i> (strain ATCC 14580 / DSM 13 / JCM 2505 / NBRC 12200 / NCIMB 9375 / NRRL NRS-1264 / Gibson 46)
Purity	≥85% (SDS-PAGE)
Sequence	MNFDSIIKNG LVILENGEQE VEVGIKDGKI AAIGQDLGTS AKIIDAKGSI VSPGMIDAHV HITEPGGGYR DEWEGYDTGT RAAAKGGVTT FIEMPLNQVP ATVDEESIQE KFKAGKGKLS VDVASYGGLV PFDLDGGIQE LDSNGVVAYK CFLATCGDRS IEGDFMNVDD YSLYEGMKQI AKTGKILSIH AENAAITDKL GEIASKNGET SLRAYVDSRP VFTEVEPIRK IILFAKETGC RVHIVHIACE EGVDEIVKAQ QEGVDITCET CTHYLYFYKE ELDHIGPVVK CSPPIREQSR LEGMWNRVLN GDISFVTS DH SPCTPDLKAT DNAFEAWGGI AGLQNNVDVL FDEGVQKRN M PLSKFAAIIA TNP AKRFNLA SKGSI AVGKD ADFVFIKKDA PYELKAEDLA YRHKISPYIG RKIGAQVVK T ILRGVEIYDK EKGISNEKVG RFIP
Source	Mammalian cell
Target Names	allB
Protein Names	Recommended name: Allantoinase EC= 3.5.2.5 Alternative name(s): Allantoin-utilizing enzyme
Expression Region	1-454
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.