



# Recombinant tRNA dimethylallyltransferase (miaA)

<b>Product Code</b>	CSB-EP714737BUI-B
<b>Abbreviation</b>	miaA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q65ZV3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Borrelia baviensis (strain ATCC BAA-2496 / DSM 23469 / PBi) (Borrelia baviensis)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MKEDRIVFIF GPTAVGKSNI LFHFPKNKAE IINVDSIQVY KEFNIASSKP SKNLMKHIKH HLVDFLDPEK EYTLGIFYEQ ALRIVKEIRQ KKKIPIFVGG TAFYFKHLKD GFPSTPLVTS EIRIYVNNLL YLKGKSYLLK ELKNVDPIRF NMLNKNDIYR IKRSLEVYYQ TGIPISQFEK KQNSEFKNIL IIGLKRSEFED LKTRISIRIN EMLNSGLLSE IKGLFSKGYN ENTPAFKGIG YNEFLWKS PCYSLNDIIS LINKNSFLYA KRQMTFFAKI SNVLWFHPED DLVSILNLIF NVDKEI
<b>Source</b>	E.coli
<b>Target Names</b>	miaA
<b>Protein Names</b>	Recommended name: tRNA dimethylallyltransferase EC= 2.5.1.75 Alternative name(s): Dimethylallyl diphosphate:tRNA dimethylallyltransferase Short name= DMAPP:tRNA dimethylallyltransferase Short name= DMATase Isopentenyl-diph
<b>Expression Region</b>	1-306
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.