



Recombinant Ribosomal RNA small subunit methyltransferase I (rsml)

Product Code	CSB-BP714754BUI
Abbreviation	rsml
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q661J3
Product Type	Recombinant Protein
Immunogen Species	Borrelia baviensis (strain ATCC BAA-2496 / DSM 23469 / PBI) (Borrelia baviensis)
Purity	≥85% (SDS-PAGE)
Sequence	MLYIVGTPIG NLEDITYRAV DVLKSVNVIF AEDTRVTIKL LSRYKINKKM ISCNVAVTENK KIRLLLDYLA KGDSVAFVSD AGTPGLSDPG SLLVAAAFRE GYKVCPIPGV SSFNTIVSVN PFRDKSVLFE GFLPNKGLKR FKKIAELYKR GDAFVLLESG HRLKLLVEI SSVSLDAKIL IGREMTKIYE EYQIGKPLEL KKYFESSKEK VKGEFTILVS RSSS
Source	Baculovirus
Target Names	rsml
Protein Names	Recommended name: Ribosomal RNA small subunit methyltransferase I EC=2.1.1.198 Alternative name(s): 16S rRNA 2'-O-ribose C1402 methyltransferase rRNA (cytidine-2'-O-)-methyltransferase Rsml
Expression Region	1-224
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.