



Recombinant Coxsackievirus B3 Genome polyprotein

Product Code	CSB-EP714758CFAW-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66282
Product Type	Recombinant Protein
Immunogen Species	Coxsackievirus B3 (strain Woodruff)
Purity	≥85% (SDS-PAGE)
Sequence	GAQVSTQKT GAHETGLNAS GNSIIHYTNI NYYKDAASNS ANRQDFTQDP SKFTEPVKDI MIKSLPALNS PTVEECGYSD RVRISITLGNS TITTQECANV VVGYGWPDY LKDSEATAED QPTQPDVATC RFYTLDSVQW QKTSPGWWWK LPDALSNLGL FGQNMQYHYL GRTGYTIHVQ CNASKFHQGC LLVVCVPEAE MGCATLNNTP SSAELLGGDS AKEFADKPVA SGSNKLVQRV VYNAGMGVGV GNLTIIPHQW INLRTNNSAT IVMPYTNSVP MDNMFRHNNV TLMVIPFVPL DYCPGSTTYV PITITIAPMC AEYNGLRLAG HQ
Source	E.coli
Protein Names	Recommended name: Genome polyprotein Cleaved into the following 12 chains: 1. Protein VP0 Alternative name(s): VP4-VP2 Protein VP4 Alternative name(s): P1A Virion protein 4 Protein VP2 Alternative name(s): P1B Virion pr
Expression Region	2-332
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.