



Recombinant 3-ketoacyl-CoA thiolase (fadI)

Product Code	CSB-MP714808YAH
Abbreviation	fadI
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q668V0
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MSKPLPLVTR QGDRIVVNG LRTPFAKQAT AYHGVPVAVDL GKIVVSELLA RSGISSELID QLVFGQVVQM PEAPNIAREI VLGTGMSVHT DAYSVSRACA TSFQAVANVA ESIIAGSVDI AIAGGADSSS VLPIGVSKAL ARTLVDANKA RSLSQKLKLF SRLRLRDLLP VAPAVAEYST GLRMGDTAEQ MAKTYGISRE DQDALALRSH QLAAEAWQQG WLHDEVMTAY IPPYREAIE DNNIRKDSTL AQYAKLRPAF DRQHGSVTAA NSTPLTDGAA AVLMMSESKA KALGLPPLGY LRSFAFSAID VWQDMLLGPS YATPLALDRA GITLADTLI DMHEAFAAQ LANLKMFAFD TFAREKLGRS QAIGEVDMSK FNVLGGSIA YGHPFAATGAR MITQTLNELR RRGGLGLTT ACAAGGLGAA MILEVE
Source	Mammalian cell
Target Names	fadI
Protein Names	Recommended name: 3-ketoacyl-CoA thiolase EC= 2.3.1.16 Alternative name(s): ACSs Acetyl-CoA acyltransferase Acyl-CoA ligase Beta-ketothiolase Fatty acid oxidation complex subunit beta
Expression Region	1-436
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.