



Recombinant Transaldolase (tal)

Product Code	CSB-BP714849YAH
Abbreviation	tal
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66ET5
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MTDKLTSLRQ ITTVVADTGD IAAMKLYQPQ DATTNPSIIL NAAQIPEYRK LIDEAIAWAR EQSSDHAQQI IDATDKLAVN IGLEILKLIP GRISTEVDAR LSYDTVASVT KAKRLIKLYN EAGISNDRIL IKLASTWQGI RAAEQLEKEG INCNLTLIFS FAQARACAEA GVFLISPFVG RILDWYKANG DQKEFAPSED PGVSVTEIY QYYKKHGYKT VVMGASFRNL GEIIELAGCD RLTIAPSLLK ELAESGPVE RKLAYTGEIQ AKPAPLTEAE FYWQHNQDPM AVDKLADGIR KFAIDQGKLE KMISDLL
Source	Baculovirus
Target Names	tal
Protein Names	Recommended name: Transaldolase EC= 2.2.1.2
Expression Region	1-317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.