



Recombinant Danio rerio Probable methylthioribulose-1-phosphate dehydratase (apip)

Product Code	CSB-MP714890DIL
Abbreviation	apip
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66I75
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSSVVNAAYA EYNGKEKGDQ EDPRVLIPQL CRLFYELGWV TGTGGGISLR HGEHIYIAPS GVQKERIQPE DLFVCDIDEK DISCPPPQKK LKKSQCTPPF MNAYTMRGAQ AVIHTHSKSA VMATLLFPGK EFRITHQEMI KGIRKGNST NFRYDDTLVV PIIENTPEEK DLKERMARAM DMYPDSCAVL VRRHGVYVWG ETWEKAKTMC ECYDYLFDIA VQMKQSGLDP SAFPTTEEKGI V
Source	Mammalian cell
Target Names	apip
Protein Names	Recommended name: Probable methylthioribulose-1-phosphate dehydratase Short name= MTRu-1-P dehydratase EC= 4.2.1.109 Alternative name(s): APAF1-interacting protein homolog
Expression Region	1-241
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.