



Recombinant Mouse Protein FAM170A (Fam170a)

Product Code	CSB-EP714908MO-B
Abbreviation	Fam170a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66LM6
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MKRRQKRKHL EIEESKEAGI SKSQEDISHP ESTGVPKAQS PGVGEVSSAS EYFSCVSSPQ KLIHRSKGTW KLLQDSSKPR SPLDQVPEGE ATTAPSQQAS SSCPSYKTCV SSLCMNKEER GMKIYYMQVQ MKKGVAISWD TKETSESLEK QPRMEEATLP EGVWVGTPPS DVSTRNLLSD SEPIGEEKEH EEKPESDSPP GSPAVEERPR AKTPDWLVTM ENGFRCMACC RVFATMESLQ EHVQYGIREG FSCHVFHLM AQLIGSMESE STQEEEEEDHT EETEKPKEEK AEEQQPTEED VGMKKPWSQC PGCVFDSPKD RRRRKDHCDN SGS
Source	E.coli
Target Names	Fam170a
Protein Names	Recommended name: Protein FAM170A Alternative name(s): Zinc finger domain-containing protein Zinc finger protein ZNFD
Expression Region	1-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.