



Recombinant *Macaca fascicularis* Mannose-binding protein C (MBL2)

Product Code	CSB-EP714911MOV-B
Abbreviation	MBL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66S50
Product Type	Recombinant Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Purity	≥85% (SDS-PAGE)
Sequence	ETVTCEDSQK ICPAVIACNS PGINGFPGKD GRDGTKGEKG EPGQGLRGLQ GPPGKLGPPG NPGSSGSPGP KGQKGDGPES PDCESSLAAS ERKALQTEMA RIKKWLTFSL GRQVGNKFFL TNGEMMTFDK VKALCAEFQA SVATPRNAE NRAIQNLIKE EAFLGITDEN TEGEFVDLTG NKLTYTNWND GEPNAGSNE DCVLLLKNGK WNDIPCSSH LALCEFPI
Source	<i>E.coli</i>
Target Names	MBL2
Protein Names	Recommended name: Mannose-binding protein C Short name= MBP-C Alternative name(s): MBP1 Mannan-binding protein Mannose-binding lectin
Expression Region	21-248
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.