



Recombinant 3-dehydroquinate synthase (aroB)

Product Code	CSB-YP714964SRW
Abbreviation	aroB
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q67N10
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MELGEERYPI YVGAGLMAQL GALVRRHLPE TRRALLVTDA NVAPLYGEAA RAALEQAGIT TARVTVPAGE ASKSLSQAYS LYSSCVRAEL DRTSAVVALG GGVVGDLAGF VAATYLRGIP LVQVPTLLA QVDSSIGGKT GVDLPQGKNL VGAFHQPSLV VADVETLKSL PPRELSAGMA EVVKHGVIIRD EEFLKYVEVQ VGNVLAGDPA VLERVVAESC RIKAEVVAAD PREKGLRAIL NFGHTVGHAL EAAIGFRWLH GECVAVGMVA AAHIARHSGI AQEGRLELRL TRLLERLDLP TTLPEGMDPK DLEPFLRRDK KIKSGVIHWV MPVRPGEVVV TPDVSPDAIR WGLEALRR
Source	Yeast
Target Names	aroB
Protein Names	Recommended name: 3-dehydroquinate synthase EC= 4.2.3.4
Expression Region	1-358
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>