



Recombinant Sulfate adenylyltransferase (sat)

Product Code	CSB-BP714978SRW
Abbreviation	sat
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67QB5
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MATQIAPHGG RLVDRWLRGP AREEALERAR RLPRVRLDAR EAADLEMIGD GALSPLTGFM GQADYRSVVA EMRLASGLLW ALPVT LAVSR AEAESIREGE EIALEDPGGR LMAVMRVAER FAYDRGAEEA RCYGT TDP AH PGVRRLLRQG EVYLGGEVWL LDRPPAPFAE YRLTPAETRA EFARRGWRTV VGFQTRNPVH RAHEYIQKCA LEICDGLLLH PLVGETKDDD LPAAVRMRAY EAILEGYFPR ERILLAVFPA AMRYAGPREA VWHALCRKNY GCTHFIVGRD HAGVGSFYGP YDAQRIFDHL DPAELGITPL FFDHTFWCRT CGAMASPKTC PHGPEARVAL SGTRVREMLY RGEAPPPEFT RPEVARVLME GLQAAARPQS VET
Source	Baculovirus
Target Names	sat
Protein Names	Recommended name: Sulfate adenylyltransferase EC= 2.7.7.4 Alternative name(s): ATP-sulfurylase Sulfate adenylate transferase Short name= SAT
Expression Region	1-393
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.