



Recombinant Putative phosphotransferase STH585 (STH585)

Product Code	CSB-MP714981SRW
Abbreviation	STH585
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67RX3
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MVAKPVIFVV SDSVGETADL VTRAALSQWG QRSVEIRRF MIDDPQSLRE VIDAARQQPT LIVYTIIVPE IRKELQRLVD MWDIPAVDIM GPMLDGISRL LGAPEKLEPG RIHRLDDDYF RRIEAVEFAV KYDDGKDPRG FLRADAVLIG VSRSSKTPVS MYLAHRRYKV ANLPLVPEVS LPKELFLVPP HKVVGLRVSP EKLHQIREER IKTIGLRSDA NYSSMERILQ ELDYAESVFK KVGCAVIDVT NKAVEETAVR VLEIINRGVN FGD
Source	Mammalian cell
Target Names	STH585
Protein Names	Recommended name: Putative phosphotransferase STH585 EC= 2.7.-.-
Expression Region	1-273
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.