



Recombinant Trigger factor (tig)

Product Code	CSB-EP714985SRW
Abbreviation	tig
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67SK1
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MKATWDKLEN NWMQFEVEVD ASEFSKSVEA AFRKMNRQVT IPGFRRGKAP RAVFERVYGK EALVQEAVDQ ILPRAYSEAL EQGQIEPISQ PEIEVVQAED GKPLIFKGKV QVTPEVTLGR LSGFGIEKQV PEVTDQEEVEE QLSALRERLA TLVTDESGEV GQGSFAVIDF EGFIDGEPFE GKGEGYTVVE IGSGTFIPGF EEGLVGAKAG ETREVTVTFP EDYHAQHLAG KQAVFKVTVH EVKKKELPEL NDEFAQQVSP FKTVEELRAD IKNRLSEAAA QKAEEDFRNK VVEAVADDAS VEVPEVLVHD KVHDMIHDFE HQLAQQGITL DMWHQITGKT HEDLHKELEP QATKLVKVDL VLA AVAKQVG LTVSDAELEA EFDRLAMYP KQQSYIRQLR ASAAYRAQLR RDLLKQKTVE HLVELNSAA
Source	E.coli
Target Names	tig
Protein Names	Recommended name: Trigger factor Short name= TF EC= 5.2.1.8 Alternative name(s): PPIase
Expression Region	1-429
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.