



# Recombinant *Oryza sativa* subsp. japonica Homeobox-leucine zipper protein HOX11 (HOX11)

<b>Product Code</b>	CSB-YP714991OFG
<b>Abbreviation</b>	HOX11
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q67UE2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MELGLSLGEA MADAGRELVL GLGMGRREEA AEAGRRDHEV RRELEFGSMS SRCGGSSPEP TVRLTLLPMV PGLGLPWPPP PPPSSSRHL EASTRGFDVN RPPSSGGGGG GGGAEQQDD VAGAALSSSP NNSAGSFPMD DFSGHGLGGN DAAPGGGGGD RSCSRASDED DGGSARKKLR LSKEQSAFLE ESFKEHSTLN PKQKLALAKQ LNLRPRQVEV WFQNRRTARTK LKQTEVDCEY LKRCCETLTE ENRRLQKELA ELRALKTVHP FYMHLPATTL SMCPSCERVA SNSAPATASS AATSSTAAPP AAPSSGGIAA TSSSAAAAAA PDHRPSSFAA LFSSPRGFPL SVAPQAQPPT SS
<b>Source</b>	Yeast
<b>Target Names</b>	HOX11
<b>Protein Names</b>	Recommended name: Homeobox-leucine zipper protein HOX11 Alternative name(s): HD-ZIP protein HOX11 Homeodomain transcription factor HOX11 OsHox11
<b>Expression Region</b>	1-362
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.