



Recombinant Human Protein PMS2CL (PMS2CL)

Product Code	CSB-MP715011HU
Abbreviation	PMS2CL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68D20
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MHAADLEKPM VEKQDQSPSL RTGEEKRDVS ISRLREAFSL RHTTENKPHS PKTPEPRRSP LGQKRGMSST STSDAISDRG VLRPQKEAVS SSQGPSDPTD RAEVEKDSGH GSTSVDSEGF SIPDTGSHCS SECVASTPGD RGSQEHVDSQ EKAPETDDSF SDVDCHSNQE DTGCKFQVLP QPTNLTSPNT KVF
Source	Mammalian cell
Target Names	PMS2CL
Protein Names	Recommended name: Protein PMS2CL Alternative name(s): PMS2-C terminal-like protein
Expression Region	1-193
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.