



Recombinant Rickettsia typhi GTPase Era (era)

Product Code	CSB-BP715093RNE
Abbreviation	era
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q68XY6
Product Type	Recombinant Protein
Immunogen Species	Rickettsia typhi (strain ATCC VR-144 / Wilmington)
Purity	≥85% (SDS-PAGE)
Sequence	<pre>MINIQKTIS VCIIGRPNSG KSTLLNRIIG EKLSIVTPKV QTTRSITGI ITLKDTQIIL YDTPGIFEPK GMLEKAMVRC AWSSLYSADL VLSIIDSLKP LDDMAHNILN QFCLLNIVPI FLLNKIDIES KYLNDIKAFL KISHPKSLLF PISALCGKNV DVLLKYIKSK AKVSPWLYAD DDITNLPMRF IAAEITREQL FLNLQQELPY KLTVQTEKFE ELKDKSIKIN QVIVISRESY KAILGKNGT KIKDIGVKSR IQMEQLFCVP VHLFLFVKVH ALWENKQEFY QYMKI</pre>
Source	Baculovirus
Target Names	era
Protein Names	Recommended name: GTPase Era
Expression Region	1-295
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>