



Recombinant *Oryza sativa* subsp. *japonica* MADS-box transcription factor 55 (MADS55)

Product Code	CSB-BP715118OFG
Abbreviation	MADS55
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q69TG5
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MARERREIRR IESAAARQVT FSKRRRGLFK KAEELAVLCD ADVALVVFSS TGKLSQFASS NMNEIIDKYT THSKNLGKTD KQPSIDLNFF LIILLRTYTN SYAYIHLLQ LEHSCSSLN EQLAEASLQL RQMRGEELEG LSVEELQQME KNLEAGLQRV LCTKDQQFMQ EISELQRKGI QLAEENMRLR DQMPQVPTAG LAVPDTENVL TEDGQSSSESV MTALNSGSSQ DNDDGSDISL KLGLP
Source	Baculovirus
Target Names	MADS55
Protein Names	Recommended name: MADS-box transcription factor 55 Alternative name(s): OsMADS55
Expression Region	1-245
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.