



# Recombinant Mouse RNA-binding protein MEX3B (Mex3b)

<b>Product Code</b>	CSB-EP715124MO
<b>Abbreviation</b>	Mex3b
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q69Z36
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MPSSLFADLE RNGTHGGGGGG GGGGGGGGSG GGETLDDQRA LQLALDQLSL LGLDSDEGAS LYDSEPRKKS VNMTECVVPV SSEHVAEIVG RQGGSGRDGD RRGFSISPTP SLEPWLPGCK IKALRAKTNT YIKTPVRGEE PVFVVTGRKE DVAMARREII SAAEHFSMIR ASRNKNTALN GAVPGPPNLP GQTTIQVRVP YRVVGLVVGP KGATIKRIQQ QTHTYIVTPS RDKEPVFEVT GMPENVDRAR EEIEAHIALR TGGIIELTDE NDFHANGTDV GFDLHHGSGG SGPGSLWSKP TPSITPTPGR KPFSSYRNDS SSSLGSASTD SYFGGGTSGS AAATSRLADY SPPSPALSFA HNGNNNNNGN GYTYTAGEAS VPSPDGGPEL QPTFDPAAP PPGTPLLWAQ FERSPGGGSA APVSSSCSSS ASSSASSSSV VFPGGGASST PSNANLGLLV HRRLHPGTSC PRLSPPLHMA TGAGEHHLAR RVRSDPGGGG LAYAAYANGL GTQLPGLPSS DTSGSSSSSS SSSSSSSSSS GLRRKGSRDC SVCFESEVIA ALVPCGHNLF CMECANRICE KSEPECPVCH TAVTQAIRIF S
<b>Source</b>	E.coli
<b>Target Names</b>	Mex3b
<b>Protein Names</b>	Recommended name: RNA-binding protein MEX3B Alternative name(s): RING finger and KH domain-containing protein 3
<b>Expression Region</b>	1-601
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.