



Recombinant Mouse RNA-binding protein MEX3B (Mex3b)

Product Code	CSB-YP715124MO
Abbreviation	Mex3b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q69Z36
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MPSSLFADLE RNGSGGGGGG GGGGGGGGSG GGETLDDQRA LQLALDQLSL LGLDSDEGAS LYDSEPRKKS VNMTECVVPV SSEHVAEIVG RQGGSGRDGD RRGFSISPTP SLEPWLPGCK IKALRAKTNT YIKTPVRGEE PVFVVTGRKE DVAMARREII SAAEHFSMIR ASRNKNTALN GAVPGPPNLP GQTTIQVRVP YRVVGLVVGP KGATIKRIQQ QTHYIVTPS RDKEPVFEVT GMPENVDRAR EEIEAHIALR TGGIILTDE NDFHANGTDV GFDLHHGSGG SGPGSLWSKP TPSITPTPGR KPFSSYRND SSSLGSASTD SYFGGGTSGS AAATSRLADY SPPSPALSFA HNGNNNNNGN GYTYTAGEAS VPSPDGGPEL QPTFDPAAP PPGTPLLWAQ FERSPGGGSA APVSSSCSSS ASSSASSSSV VFPGGGASST PSNANLGLLV HRRHLPGTSC PRLSPPLHMA TGAGEHHLAR RVRSDPGGGG LAYAAAYANGL GTQLPGLPSS DTSGSSSSSS SSSSSSSSSS GLRRKGSRDC SVCFESEVIA ALVPCGHNLF CMECANRICE KSEPECVPCH TAVTQAIRIF S
Source	Yeast
Target Names	Mex3b
Protein Names	Recommended name: RNA-binding protein MEX3B Alternative name(s): RING finger and KH domain-containing protein 3
Expression Region	1-601
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.