



Recombinant Propionibacterium acnes Elongation factor Tu (tuf)

Product Code	CSB-MP715142PSG
Abbreviation	tuf
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6A6L7
Product Type	Recombinant Protein
Immunogen Species	Propionibacterium acnes (strain KPA171202 / DSM 16379)
Purity	≥85% (SDS-PAGE)
Sequence	MAKAKFERTK PHCNIGTIGH IDHGKTTTLTA AISKVLHDKY PELNEESPFQ QIDKAPEERQ RGITISIAHI EYQTEKRHYA HVDCPGHADY VKNMITGAAQ MDGAILVVAA TDGPMPQTRE HVLLARQVGV PAIVVALNKC DMVDDEELIE LVEMEVRELL TSQEFDGDNC PVVRISAFQA LQGDEKWTQS ILDLMDAVDE YIPQPERDLQ KPFLMPIEDV FTITGRGTVV TGRVERGVVK TGEEVEIVGI HEKTQKTTVT GVEMFRKILD EGRAGENVGV LLRGTKKEDV VRGMVLSKPG STTPHTDFEG QVYVLKKDEG GRHKPFFSHY SPQFYFRTTD VTGTVELPEG TEMVMPGDNT DMTVHLIHPV AMEDQLKFAI REGGRTVGAG RVTKIHK
Source	Mammalian cell
Target Names	tuf
Protein Names	Recommended name: Elongation factor Tu Short name= EF-Tu
Expression Region	1-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.