



# Recombinant Leifsonia xyli subsp. xyli Exodeoxyribonuclease 7 large subunit (xseA)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP715184LAAC-B                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Abbreviation</b>      | xseA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                               |
| <b>Uniprot No.</b>       | Q6ADU9                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Immunogen Species</b> | Leifsonia xyli subsp. xyli (strain CTCB07)                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Sequence</b>          | MRGTRVTETA SAPRMAPGPP TLDDPWPVAL LASKIRGWIE RLGTAWVEGE<br>ITQWGVSGGN VYGKLDLNE DATVGFTIWS SVKARIPADL KQGDRVIAAV<br>KPNYWLKGGT LTMQVSDMRH VGLGDLLERL ERLRAQLRAE GLFRPERKKR<br>LPFLPHTIGL VTGKDSDAEK DVLRNAQLRW PQVRFRTVYA AVQGDRTVPE<br>VTAALRELD APEVEVIIVA RGGGDFQNL GFSDSLLRA AAGLSTPLVS<br>AIGHEADRPL LDEVADLRAS TPTDAAKRVV PDVAEELVRV HQARARIGTR<br>LTHIRHEID RIGHLRTRPA LASGSWIVDS RAQDLTRFVA RGAELVERCV<br>DREARVAEL RGQLRALSPQ ATLERGYAIV QNAAGRVA A PGEAPAGTEL<br>RITVSGGSLA ATAGKALPSP AQGATNGSLA APRGK |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Target Names</b>      | xseA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Protein Names</b>     | Recommended name: Exodeoxyribonuclease 7 large subunit EC= 3.1.11.6<br>Alternative name(s): Exodeoxyribonuclease VII large subunit Short name= Exonuclease VII large subunit                                                                                                                                                                                                                                                                                                                                     |
| <b>Expression Region</b> | 1-435                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                      |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life                                                                                                                                                                                                                                                                                |



of lyophilized form is 12 months at -20°C/-80°C.