



Recombinant Danio rerio Protein FAM199X (fam199x)

Product Code	CSB-MP715254DIL
Abbreviation	fam199x
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AXJ7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MSEGVYEKFL APEEPFPLLS QRGNFSEDAT LDVSDFGCQL SSCHRTDPLH RFHSSRWNLT SCGTSVASSE CSEELFSSVS VGDQDDCYSL LDDQELTSFD LFPEGSVCSD VSSSISTYWD WSDSEFEWQL PGSDIASGSD VLSDIIPSVP SSPCLVSKRK SKPHRNLDL PWSAMTNDEQ VEYIEYLSRK VSTEMGLREQ LDIIKIDPS AQILPTDSEF IELNCLTDE KLKQVRSYIR EHSPRQRPN TRDNWKRSM SNGSTSGVSA HSSSNASMVS STSSSTASTG SNSSTNISRA HSDGNLATAA ERIRDSKKRS KQRKMQQKAL RKRQLKEQRQ ARKERLSGLF MNEEVLSLKV TEEEDHCDDV DVLM
Source	Mammalian cell
Target Names	fam199x
Protein Names	Recommended name: Protein FAM199X
Expression Region	1-374
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.